

Monitoring of Cold Chain and Production

DATA COLLECTION ACCORDING TO GXP AND FDA GUIDELINES



11471 台北市内湖區新明路273巷6號1樓 1F., No.6, Ln. 273, Xinming Rd., Nelhu Dist., Talpel City 11471, Talwan (R.O.C.) 服務專線Tel : (02)8792-3722

服務傳真Fax: (02)8792-3761

電子信箱Email: info@grandever-biotech.com.tw 公司網址Website: www.grandever-biotech.com.tw







To ensure your quality

By using the right measuring devices, a high and constant product quality can be guaranteed. This often requires checks in the incoming goods department, during processing, on the finished end product as well as during storage and logistics. Xylem Analytics, offers you the necessary measuring devices that are suitable for use in the GxP environment, as well as software solutions for your diverse applications. You can find a selection of applications and devices here:

Which device is right for you?

| | TMX 3x0/4x0 | EBI 20 Series |
|---|-------------|---------------|
| Measurement | | |
| For direct measurement | 1 | 1 |
| For monitoring | | ✓ |
| Cool Chain Control | | |
| Refrigerator temperature control | ✓ | ✓ |
| Warehouse monitoring | | ✓ |
| Logistics | | ✓ |
| Mappings | | ✓ |
| Laboratory | | |
| Clinical Trails | ✓ | ✓ |
| Temperature control in refrigerators | ✓ | ✓ |
| Incubators | ✓ | ✓ |
| Process Control | | |
| Process validation | | |
| Qualification | | |
| Steamsterilization | | |
| H ₂ O ₂ and EtO Sterilization | | |
| Depyrogenation / heat tunnel | | |



For more information in our catalogues.

Visit us at **www.ebro.com** for many more products!



| EBI 300 Series | EBI 25 Series | EBI 310 Series | EBI 11 Series | EBI 12 Series |
|----------------|---------------|----------------|---------------|---------------|
| | | | | |
| ✓ | | ✓ | | |
| ✓ | ✓ | ✓ | ✓ | ✓ |
| | | | | |
| ✓ | ✓ | ✓ | | |
| ✓ | ✓ | ✓ | | |
| ✓ | | ✓ | | |
| ✓ | ✓ | ✓ | ✓ | ✓ |
| | | | | |
| ✓ | ✓ | ✓ | ✓ | ✓ |
| ✓ | ✓ | ✓ | | |
| ✓ | ✓ | ✓ | | |
| | | | | |
| | | | ✓ | ✓ |
| | | ✓ | ✓ | ✓ |
| | | | ✓ | ✓ |
| | | | | ✓ |
| | | | | ✓ |
| | | | | |





Data Logger to measure during manufacturing

EBI 11 Series



Mappings, clinical trails, process validation and qualification and steam sterilization

The EBI 11 offers:

- Small data logger for applications with little space
- 16.5 mm diameter
- Measuring range of -30 °C to +150 °C
- Logger to measure temperature and / or pressure available
- Calibration certificate included
- Stainless steel (V4A) housing
- Sets with interface, software and accessories for various applications are available.
- Winlog.pro software for evaluation
- Waterproof housing (IP68)
- Prepared IQ / OQ documents are available for Winlog.pro software
- FDA 21 CFR part 11 compliant



大久生物科技股份有限公司 GRANDEVER BIOTECHNOLOGY CO., LTD.

Process data logger

EBI 12 Series



Mappings, clinical trails, process validation and qualification, steam sterilization, H_2O_2 and EtO sterilization and depyrogenation / heat tunnel

The EBI 12 offers:

- Logger to measure temperature, pressure, humidity, conductivity and precise vacuum
- Large selection of different types
- Calibration certificate included
- Stainless steel (V4A) housing
- Sets with interface, software and accessories for various applications are available.
- TÜV certified Winlog.validation software for evaluation
- Waterproof housing (IP68)
- Prepared IQ / OQ documents are included in the Winlog.validation software package
- FDA 21 CFR part 11 compliant





入久生物科技股份有限公司 GRANDEVER BIOTECHNOLOGY CO., LTD.

Digital Refrigerator Thermometer

TMX 3x0 and 4x0 Min/Max Alarm Thermometer with external probe



For temperature control in refrigerators, laboratory and clinical trails.

The TMX 3x0 / 4x0 offers:

- Observe two temperature zones with one device
- Easy usage
- Reliable control
- Acoustic alarm for active alarm management
- Optical alarm to see how long the problem occurs





Standard Cold Chain Data Logger

EBI 20 Series Robust Basic Data Logger



For temperature control in refrigerators, warehouse monitoring and logistic control, mappings, laboratory and clinical trails.

The EBI 20 offers:

- Data logger versions for temperature and humidity measurements available
- With external or internal temperature probe
- LED shows if a limit has been exceeded
- Measurement range of -30 °C to +100 °C
- User replaceable battery
- 2 point calibration certificate included
- Interface and Desktop software are necessary
- Winlog.basic software free of charge
- Waterproof housing (IP 67) except for humidity logger



STANDEVER BIOTECHNOLOGY CO., LTD.

PDF - Cold Chain Data Loggers with USB connector

USB Data Logger EBI 3x0 Series Continuous temperature monitoring sensitive goods



For temperature control in refrigerators, warehouse monitoring and logistic control, mappings, laboratory, clinical trails and for device qualification.

The EBI 3x0 offers:

- USB-A connector enables you to read the logger out on every windows based PC worldwide without the need of an interface
- Automatic PDF generation without software
- Configuration via the website <u>www.ebi300.com</u> or the free desktop software Winlog.basic
- LED shows if a limit has been exceeded
- Measurement range of -30 °C to +70 °C
- Several external penetration probes available
- Digital humidity probes available
- User replaceable battery
- 3 point calibration certificate included
- Waterproof housing (IP 65)



大久生物科技股份有限公司 GRANDEVER BIOTECHNOLOGY CO., LTD.

Wireless Cold Chain Monitoring System

EBI 25 Series Automated data collection and processing



For temperature control in refrigerators, warehouse monitoring, climatic chambers and stability spaces, laboratory, incubators, mappings and clinical trails.

The EBI 25 offers:

- Data logger versions for temperature and humidity measurements available
- With external or internal temperature probe
- Measurement range of -200 °C to +199.9 °C
- User replaceable battery
- Alert notification by email
- 2 point calibration certificate included
- Interface and webbased software are necessary
- Winlog.web software
- FDA 21 CFR part 11 compliant
- Waterproof housing (IP 67) except for humidity logger



KompetenzCentrum ebro®

Theory and practice are the key to our successful transfer. Our instructors are all experts in their own fields.

The aim is to achieve an in-depth understanding of hardware and software alike.

www.ebro.com/en/training

Founded in 1968, in its more than 50-year history, ebro® has specialized in high-quality measurement technology such as handheld instruments, thermometers and data loggers. Are you looking for devices and software to measure temperatures, humidity, pH, oil qualities, pressure and Brix? Then you have found your partner with ebro®. From measurement to recording to evaluation, we offer you quality products and solutions tailored to your requirements.





11471 台北市内湖區新明路273巷6號1樓 1F., No.6, Ln. 273, Xinming Rd., Nelhu Dist, Taipel City 11471, Taiwan (R.O.C.) 服務專線Tel: (02)8792-3722

服務專線Tel: (02)8792-3722 服務傳真Fax: (02)8792-3761

wsh博具Fax:(02)8/92-3/61 電子信箱Email:info@grandever-biotech.com.tw 公司網址Website: www.grandever-biotech.com.tw





EBI 25 Wireless System

MONITORING OF TEMPERATURE, HUMIDITY AND MANY OTHER MEASUREMENTS





11471 台北市内湖區新明路273巷6號1樓 1F., No.6, Ln. 273, Xinming Rd., Neihu Dist., Taipel City 11471, Taiwan (R.O.C.) 服務專線Tel: (02)8792-3722

服務傳真Fax: (02)8792-3761

電子信箱Email: info@grandever-biotech.com.tw 公司網址Website: www.grandever-biotech.com.tw







Online Monitoring of Temperature, Humidity and Many Other Measurements

The EBI 25 system for wireless monitoring of temperature, humidity and other measurements assures that perishable goods are produced and stored at the right conditions at all times. Other measurements can be integrated using Mod bus over IP.



Continuous monitoring

of all areas where food, pharmaceutical and medical products are produced and stored.



with the continuous monitoring of all crucial temperature points – 24 hours on 365 days a year.



Quick intervention before it is too late

with a very wide and flexible alarm management (notification via email or acoustical).



Automatic documentation

through automatic recording, storage, and reporting of the measured data. Problems can easily be traced and audit reports can be generated instantly.



Worldwide access to the measurement data

via internet – at any time, any place and any PC with internet and network access. Security is ensured using a login and password.

Easy handling

saves time, effort and energy.



Applications

Food sector

- Continuous monitoring of cooling equipment, cool rooms and deep-freeze rooms
- Monitoring of room temperatures and humidity values in a store

Pharmaceutical / medical sector

- Continuous monitoring of temperature sensitive goods like blood plasma or medicine
- Temperature and humidity monitoring and other measurements in cooling chambers, clean rooms and incubators

Benefits

Food sector

- Compliance with the HACCP guidelines
- All devices according to DIN EN 12830
- Up to date reports always ready for food inspectors

Pharmaceutical / medical sector

- FDA 21 CFR Part 11 functionality
- Up to date reports generated on-demand for audits and inspections



Technical Data

| Туре | EBI 25-T | EBI 25-TE | EBI 25-TX | EBI 25-TH |
|--------------------------------|---|------------------|--|---|
| Probes | Internal | External | External for very deep temperatures | Internal: temperature External: humidity |
| Measurement range: temperature | -30 °C +60 °C | -40 °C +85 °C | -200 °C +199.9 °C | -30 °C +60 °C |
| Measurement range: humidity | | | | 0% rH 100% rH |
| Accuracy: temperature | ± 0.5 °C (-20 °C +40 °C) ± 0.8 °C for the remaining measurement range | | ±2.0 °C (-200 °C100 °C) ±0.5 °C (-100 °C20 °C) ±0.2 °C (-20 °C +60 °C) ±0.5 °C (+60 °C +199.9 °C) | ±0.5 °C (-20 °C +40 °C) ±0.8 °C for the remaining measurement range |
| Accuracy: humidity | | | | ±3 % |
| Resolution: temperature | 0.1 °C | 0.1 °C | 0.1 °C | 0.1 °C |
| Resolution: humidity | | | | 0.1% rH |
| Memory capacity | 288 measurements | 288 measurements | 288 measurements | 576 measurements (288 measurements per channel) |
| Sensor | NTC | NTC | Pt 1000 | NTC |
| Protection class | IP 67 | IP 67 | IP 67 | IP 20 |

General technical data

| Measurement interval | 1 min. to 24 hours, adjustable |
|----------------------|---|
| Radio frequency | 868 MHz |
| Battery | 3.6 V lithium (user replaceable) |
| Battery lifetime | Up to 2 years, depending on transmission rate |
| Storage temperature | -40 °C +70 °C |



USB cable / Ethernet

Functionality

The innovative EBI 25 system monitors wirelessly, transmits the data in real time to an interface (base station) and from there to the desired server or PC.



EBI 25 series data loggers

- Precise measurement of temperature and humidity (depends on logger type)
- Very large range of up to 500 m in a free field
- Long battery lifetime
- Easy installation

Base station: IF 400 interface

- Collects and stores the data of all connected EBI 25 data loggers
- Connection of up to 50 loggers per interface possible
- Stores up to 576 measurements per logger
- Direct connection of any number of interfaces to a PC or the network
- Audible alarm (with optional alarm box)

Embedding of other measurements possible:

• Pressure

• CO₂

• NH₃

• O₂

• Many others



For more information please contact us:

+49 841 95478-0

As soon as a temperature, humidity or any other measurement has exceeded a user defined limit, an alarm is immediately sent via email.

Evaluation on the PC









Work place

Measurement management software: Winlog.web and Winlog.wave

Print

Email distribution

Winlog.wave:

Basic version for local PC usage


- Single PC solution: The basic version for easy measurement data evaluation on a single PC- no network required.
- Flexible alarm management: graphic and email alarm notifications upon user defined conditions
- Connection to the IF 400 via USB
- FDA 21 CFR Part 11 data security functionality

Winlog.web:

Professional version internet and local network based use

- Web based client/server solution: the measurement data can be evaluated on all PCs and smartphones via the internet or connected to the local network
- Very flexible and wide alarm management: alarm notifications upon user defined conditions, alarm notification via email; visual and audible alarm via the graphical user interface
- Connection of the interface IF 400 via USB and Ethernet
- FDA 21 CFR Part 11 data security functionality
- Management of larger data sets
- IQ / OQ documentation available

Quick graphical overview of all measurement points including an image or floor plan view





Types of filters

There are two different filters for the EBI 25-TH:



The PTFE-filter AH 100 provides little protection from mechanical strain but seals the probe from dust – protection class IP 60. In dusty environments, unprotected humidity probes, or probes protected by sintered filters, can quickly be covered by dust. This massively affects the humidity measurement, or prohibits it entirely. It won't happen with the PTFE filter.



The stainless steel sintered filter AH 300 increases the protection class to IP 40. The stainless steel massively increases the mechanical protection, e.g. from crushing the filter. Nothing will happen to your probe! In addition, the stainless steel is highly resistant to corrosion even in high temperature applications.



Accessories for the EBI 25 system



AL 250 – Protection Box for EBI 25 TE and TX



EBI 2 AB-2 – Alarmbox to connect to interface IF 400



AL 251 – orange Flash Light and buzzer combination



EBI IF 400 – Interface including antenna



AG 152 – Wall mount for EBI 25



AL 116 – external antenna for the interface EBI IF 400 interface for more power, optional



How to order:

| Туре | Description | Part No. |
|------------|--|-----------|
| EBI 25-T | Wireless temperature logger (with internal sensor) | 1340-6200 |
| EBI 25-TE | Wireless temperature logger (with external probe) | 1340-6201 |
| EBI 25-TX | Temperature data logger (with probe TPX 25-3) | 1340-6204 |
| EBI 25-TX | Temperature data logger (without probe) | 1340-0025 |
| TPX 25-3 | Pt 1000 probe for EBI 25-TX, 3 m | 1341-0025 |
| TPX 25-5 | Pt 1000 probe for EBI 25-TX, 5 m | 1341-0026 |
| TPX 25-7,5 | Pt 1000 probe for EBI 25-TX, 7,5 m | 1341-0027 |
| TPX 25-10 | Pt 1000 probe for EBI 25-TX, 10 m | 1341-0028 |
| EBI 25-TH | Wireless temperature / humidity logger | 1340-6202 |
| AH 100 | PTFE filter for EBI 25-TH | 1340-5627 |
| AH 300 | Stainless steel sintered filter for EBI 25-TH | 1340-5625 |
| AL 250 | Protection Box for EBI 25 TE and TX | 1248-0250 |
| EBI 2 AB-2 | Alarm box to connect with interface IF 400 | 1613-1301 |
| AL 251 | Orange Flash Light and Buzzer combination | 1340-6233 |
| EBI IF 400 | Interface incl. antenna | 1340-6210 |
| AG 152 | Wall mount for EBI 25 | 1340-6215 |
| AL 116 | External Antenna | 1340-6211 |
| Winlog.web | Web based Evaliation Software – Server solution | 1340-2390 |









養品費目 機合問發製作,電路和自動控制配電盤之設 FDA/USP/EPA美國藥典完全無菌等級減 器—減菌確效&冷鏈倉備溫渥度給 計及製作,HMI/PIC 程式設計,系統規 菌清潔產品:異丙醇/變性乙醇。過乙酸 測盤控。精密型/工業型電子天 劃,SCADA圖控監控軟體規劃整合設計。 過氣化氫級能子劑/預稀釋過氧化氫溶 液/酚醛減菌劑/四級胺減菌劑/異丙醇 除溫機/負離子產生機。初放/中 機/酚醛減菌劑/四級胺減菌劑/異丙醇 除溫機/負離子產生機。初放/中 經/於於於於一/經/向性/經/位去/經// 到,SANA國程並在教院的到金含級引。 國被電機元件: INV終餘牌變頻器,POD人機 於酚醛滅菌劑/四級胺滅菌劑/異两醇 達、MS電磁開關、荷重元、SENSOR酸鹼度/ 樂電度/溶氣量/液(氣)體流量/壓力計。

美國原裝進口專業製藥工業ISO Class5 德國無線(有線)溫濕度壓力紀錄 效/高效 (HEPA) /化學/耐高溫濾網。基因檢測。組織染色試劑

11471 台北市内湖區新明路273巷6號1樓 1F., No.6, Ln. 273, Xinming Rd., Neihu Dist., Taipei City 11471, Taiwan (R.O.C.)

服務專線Tel: (02)8792-3722 服務傳真Fax: (02)8792-3761

電子信箱Email: info@grandever-biotech.com.tw 公司網址Website: www.grandever-biotech.com.tw





Xylem Analytics Germany Sales GmbH & Co. KG, ebro Peringerstr. 10 85055 Ingolstadt, Germany Tel +49 841 95478-0 Fax +49 841 95478-80 ebro@xyleminc.com www.xylemanalytics.com

ebro® is a registered trademark of Xylem Inc. or one of its subsidiaries. Technical changes reserved. © 2015 Xylem Analy tics Germany Sales GmbH & Co. KG.



EBI 300/310/330

USB DATA LOGGERS FOR COLD CHAIN MONITORING





Cold Chain Monitoring

The easy to use data loggers with USB connection monitor the temperature and/or humidity during transport and storage of sensitive goods like medicine, food, serums etc.. Measurement reports are created automatically as PDF files when you connect the logger to a PC.

The EBI 300 and EBI 310 USB data loggers are suitable for multi-use, the EBI 330 data loggers are single-use versions which can be ordered preconfigurated and are used especially when returning a more expensive multi-use logger to the sender after a shipment is difficult. Please contact us for more information.



Program | Measure

- Programming of the logger with the help of the free online configurator at www.ebi300.com or optional via the software Winlog.basic or Winlog.pro
- Set optional limits and start to record the measurement data



Connect | Readout

- · Connection of the logger to any PC via the USB port
- Automatic generation of a PDF report with all important measurement data



Evaluate | Archive

- · Store, save or email the PDF report
- Further processing of the measurement data with the software Winlog.basic or Winlog.pro





Benefits

- Direct USB connection
- · Automatic PDF report generation
- Programmable at www.ebi300.com no special software for programming and readout required but available
- Indication of alarm status via flashing I FD
- Data integrity
- Conforms with FDA 21 CFR Part 11 and ATP
- The data loggers help you to comply with GMP and VO (EG) 37/2005
- Free firmware updates at your place via software
- Excellent value-for-money







Digital interface

- Digital interface between loggers and external probes (at EBI 300 TE, EBI 300 TH, EBI 310 TE, EBI 310 TH and EBI 310 TX).
- Data logger functions as data collector with optional internal sensor
- Easy exchange of the external probes e.g. for calibration: remove and send probe, connect replacement probe, measure!
- No calibration of the data collector required, if internal probe is not used!



Which data logger is right for you?

Every EBI 3x0 USB data logger have the afore mentioned properties. Depending on the application, claim and your purse, there are different requirements for which we have the right devices. The following overview shall help making the decision.

| | EBI 330-T30 | EBI 300 | EBI 300 TE | EBI 300 TH | EBI 310 | EBI 310 TE | EBI 310 TX | EBI 310 TH |
|--|-------------|---------|------------|------------|----------|------------|------------|------------|
| Applications | | | | | | | | |
| Monitoring of deep temperatures | | | | | | √ | / * | |
| Monitoring of high temperatures | | | | | | 1 | / * | |
| Humidity monitoring | | | | 1 | | | | 1 |
| Storage monitoring | | 1 | 1 | 1 | 1 | 1 | / * | 1 |
| Transport monitoring | / | 1 | / | 1 | 1 | 1 | / * | 1 |
| Process monitoring | | | / | 1 | | 1 | √ * | 1 |
| Usage within dry ice | | | | | | | | |
| Measurement channels | | | | | | | | |
| Internal temperature channel | ✓ | 1 | √ | 1 | 1 | 1 | / | 1 |
| External temperature channel | | | 1 | 1 | | 1 | 2 * | 1 |
| Sensor cable | | | 1 | | | 1 | / * | |
| High precision (Pt 1000) | | | | | 1 | 1 | / | 1 |
| Humidity channel | | | | 1 | | | | 1 |
| Usage | | | | | | | | |
| Multi-use | | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Single-use | ✓ | | | | | | | |
| Calibration certificate | | | | | | | | |
| Including factory calibration certificate | | 1 | 1 | 1 | 1 | 1 | / * | / |
| Batch calibration certificate available upon request | ✓ | | | | | | | |
| Other features | | | | | | | | |
| Display | | 1 | 1 | 1 | 1 | 1 | 1 | 1 |
| Very flexible alarms (5 limits and MKT) | | | | | 1 | 1 | / | 1 |
| High memory capacity (120.000 measurements) | | | | | √ | √ | / | / |

^{*} with connected, exchangeable sensors



EBI 330 Single-Use USB Data Loggers

General technical specifications: valid for both EBI 330 data logger types*

| Memory capacity | 20,000 measurements |
|---|--|
| Alarm | 2 limits |
| PDF creation | PDF/A 1b |
| LED | Yes (red and green) |
| Resolution | 0.1 ℃ |
| Storage temperature | -40 °C +85 °C (-40 °F +185 °F) |
| Sample rate | 1 min. to 24 hours |
| Measurement modes | Start with key pressStart immediately until end of memory |
| Maximum start delay | 24 h |
| Dimensions (L x W x H) | 80 x 28 x 12 mm |
| Housing material | ABS |
| Protection class | IP 65 |
| Certificate | Batch calibration certificate available on request |
| * Please find the exact technical data of e | ach EBI 330 data logger type on this page. |
| | |

EBI 330-T30 Single-Use USB Data Logger Standard version



Technical Data Measurement rang / operating temperature -30 °C ... +60 °C (-22 °F ... +140 °F) Accuracy ±0.5 °C (-20 °C ... +40 °C / -4 °F ... +104 °F) ±0.8 °C for the remaining measurement range Sensor NTC Battery Lithium button cell (CR2032), 3V Battery lifetime 100 days







EBI 300 Multi-Use USB Data Loggers

General technical specifications: valid for EBI 300 data logger types*

| Memory capacity | 40,000 measurements | |
|---------------------|--|---|
| Alarm | 2 limits | |
| PDF creation | PDF | |
| LED | Yes (red) | |
| Storage temperature | -40 °C +85 °C (-40 °F +185 °F) | |
| Sample rate | 1 min 24 h | |
| Measurement modes | Endless measurement Start / Stop Measurement until end of memory Start with key press | |
| Display | Value, MIN / MAX, until end of memory, alarm on / off | |
| Maximum start delay | 72 h | |
| Housing material | Polycarbonate | |
| Certificate | Factory calibration certificate | |
| Battery | Lithium button cell (CR 2450), 3 V | • |
| Battery lifetime | Up to 2 years, depending on applications | |

^{*} Please find the exact technical data of each EBI 300 data logger type on this double page.

EBI 300 Multi-Use USB Data Logger

Standard version



Technical Data

| Measurement range/operating | -30 °C +70 °C (-22 °F +158 °F) |
|-----------------------------|--|
| temperature | By connecting an external probe, the temperature |
| | measurement range can be extended. |
| Accuracy | ± 0.5 °C (-20 °C +40 °C / -4 °F +104 °F) |
| | ±0.8 °C for the remaining measurement range |
| Sensor | NTC |
| Resolution | 0.1 °C |
| Dimensions (L x W x H) | 80 x 33 x 14 mm |
| Protection class | IP 65 |



The device has been certified together with the EB 4401 food inspection case.

Accessories for EBI 300, EBI 300 TE and EBI 300 TH



EBI 300 TE + EBI 300 WM2



EBI 300-WM2 Wall Mount for EBI 300 / 310



EBI 300 WM3 transportation mount for EBI 300/310 made of stainless steel



EBI 300 TE Multi-Use USB Data Logger with external temperature probe Fast, flexible core temperature measurements

To alonical Data

Dimensions (L x W x H)

Protection class



| lechnical Data | |
|--|--|
| Measurement range external temperature | -35 °C +70 °C (-31 °F +158 °F) |
| Measurement range internal temperature / operating temperature | -30 °C +70 °C (-22 °F +158 °F) |
| Accuracy (internal and external) | \pm 0.5 °C (-20 °C +40 °C / -4 °F +104 °F) \pm 0.8 °C for the remaining measurement range |
| Probe | NTC, Stainless steel, Ø 4mm, L = 50mm, pointed |
| Cable | PVC, $L = 1m$, waterproof, oilproof and food safe |
| Resolution | 0.1 °C |

91 x 33 x 14mm

- Simultaneous measurement of core temperature and ambient temperature
- Internal temperature probe usable additionally

EBI 300 TH Multi-Use USB Data Logger with external humidity and temperature probe

Technical Data

Measurement range temperature /

Relative humidity monitoring in storages and during transports



| operating temperature | |
|----------------------------|---|
| Accuracy (internal) | \pm 0.5 °C (-20 °C +40 °C / -4 °F +104 °F) \pm 0.8 °C for the remaining measurement range |
| Accuracy (external) | \pm 0.5 °C (-20 °C +40 °C / -4 °F +104 °F) \pm 1.0 °C for the remaining measurement range |
| Probe temperature | 0% rH 100% rH |
| Measurement range humidity | \pm 3 % between 10 % rH 90 % rH (at 25 °C/+77 °F) \pm 5 % for the remaining measurement range |
| Accuracy humidity | \pm 3 % between 10 %rH 90 %rH (at 25 °C/+77 °F) \pm 5 % for the remaining measurement range |
| Probe humidity | Capacitive |
| Resolution temperature | 0.1 ℃ |
| Resolution humidity | 0.1 % rH |
| Dimensions (L x W x H) | 129 x 33 x 14mm |
| Protection class | IP 20 |
| | |

-30 °C ... +70 °C (-22 °F ... +158 °F)

 Internal temperature probe usable additionally



EBI 310 Multi-Use USB Data Loggers

General technical specifications: valid for all EBI 310 data logger types*

| Memory capacity | 120,000 measurements | |
|---------------------|--|--|
| Alarm | 5 ranges | |
| PDF creation | PDF/A 1b | |
| LED | Yes (red and yellow) | |
| Storage temperature | -40 °C +85 °C (-40 °F +185 °F) | |
| Sample rate | 1 s 24 h | |
| Measurement modes | Endless measurement Start / Stop Measurement until end of memory Start with key press | |
| Display | Value, MIN / MAX, until end of memory, alarm on / off | |
| Maximum start delay | 72 h | |
| Housing material | Polycarbonate | |
| Certificate | Factory calibration certificate | |

^{*} Please find the exact technical data of each EBI 310 data logger type on on the next pages.

EBI 310 Multi-Use USB Data Logger

High precision version



Technical Data

| Measurement range/operating | -30 °C +75 °C (-22 °F +167 °F) | |
|-----------------------------|--|--|
| temperature | By connecting an external probe, the temperature | |
| | measurement range can be extended. | |
| Accuracy | ± 0.2 °C (-30 °C +30 °C / -22 °F +86 °F) | |
| | ±0.5 °C for the remaining measurement range | |
| Sensor | PT 1000 | |
| Resolution | 0.1 °C | |
| Dimensions (L x W x H) | 80 x 33 x 14 mm | |
| Protection class | IP 65 | |
| Battery | Lithium button cell (CR 2450), 3 V | |
| Battery lifetime | Up to 2 years, depending on applications | |

Accessories for the EBI 310, EBI 310 TE, EBI 310 TX and EBI 310 TH



EBI 300-WM2 Wall Mount for EBI 300 / 310



EBI 300 WM3 transportation mount for EBI 300/310 made of stainless steel



EBI 310 TE Multi-Use USB Data Logger with external precision temperature

probe

Measurement of high and low temperatures



- Simultaneous measurement of core temperature and ambient temperature
- Internal temperature probe usable additionally

| Technical Data | | |
|--|---|--|
| Measurement range external temperature | -200°C +250°C (-328 °F +482 °F) | |
| Measurement range internal temperature / operating temperature | -30 °C +75°C (-22 °F +167 °F) | |
| Accuracy (internal) | \pm 0.2 °C (- 30 °C + 30 °C / -22 °F +86 °F) \pm 0.5 °C for the remaining measurement range | |
| Accuracy (external) | ± 2.0 °C (- 200 °C100 °C / -328 °F148 °F) ± 1.0 °C (-100 °C20 °C / -148 °F4 °F) ± 0.2 °C (- 20 °C +60 °C / -4 °F +160 °F) ± 0.5 °C (+60 °C + 250 °C / +160 °F +482 °F) | |
| Probe | Pt 1000, Stainless steel, Ø 5mm, L = 50mm, blunt | |
| Cable | PTFE, L = 1m, waterproof, oilproof and food safe | |
| Resolution | 0.1 °C | |
| Dimensions (L x W x H) | 91 x 33 x 14 mm | |
| Protection class | IP 65 | |
| Battery | Lithium button cell (CR 2450), 3V | |
| Battery life time | Up to 2 years, depending on applications | |
| | | |





EBI 310 TH Multi-Use USB Data Logger with with external humidity and temperatore probe

Relative humidity monitoring in storages and during transports



| -30 °C +75°C (-22 °F +167 °F) |
|---|
| $\pm~0.2~^{\circ}\text{C}$ (- $30~^{\circ}\text{C}$ \dots + $30~^{\circ}\text{C}$ / - $22~^{\circ}\text{F}$ \dots +86 $^{\circ}\text{F}$) $\pm~0.5~^{\circ}\text{C}$ for the remaining measurement range |
| \pm 0.5 °C (0 °C + 60 °C / +32 °F +140 °F) \pm 0.8 °C for the remaining measurement range |
| Pt 1000 |
| 0 % rH 100 % rH |
| \pm 2 %between 10% rH 90 %rH (at 25 °C /+77 °F) \pm 4 % for the remaining measurement range |
| capacitive |
| 0.1 °C |
| 0.1 % rH |
| 129 x 33 x 14mm |
| IP 20 |
| Lithium button cell (CR 2450), 3V |
| Up to 2 years, depending on applications |
| |

 Internal temperature probe usable additionally

EBI 310 TX Multi-Use USB Data Logger with temperature-two-channel-adapter Temperature monitoring in storages and during transport, process monitoring





| Technical Data | |
|--|--|
| Measurement range external temperature | -200°C + 400°C (-328 °F +752 °F), dependent on probe type |
| Measurement range internal temperature / operating temperature | -30 °C +75°C (-22 °F +167 °F) |
| Accuracy (internal) | \pm 0.2 °C (- 30 °C + 30 °C / -22 °F +86 °F) \pm 0.5 °C for the remaining measurement range |
| Probe | Pt 1000 |
| Resolution | 0.1 °C |
| Dimensions (L x W x H) | 111 x 33 x 14mm |
| Protection class | IP 65 |
| Battery | Lithium button cell (CR 2450), 3V |
| Battery life time | Up to 2 years, depending on applications |

exchangeable sensors

- Up to two exchangeable probes can be connected; not included, see the following page
- Internal temperature probe usable additionally



Exchangeable probes for EBI 310 TX



TPX 310-P1

- Measurement range: -200 °C ... +200 °C (-328 °F ... +392 °F)
- Needle: L = 45 mm, $\emptyset = 5 \text{ mm}$, blunt
- Cable: PTFE, L = 3 m

| Temperature | | Accuracy |
|-------------|-------------|----------|
| -200100 °C | -328148 °F | 1.7 °C |
| -10020 °C | -1484 °F | 1.2 °C |
| -20+60 °C | -4+140 °F | 1.0 °C |
| +60+200 °C | +140+392 °F | 1.7 °C |



TPX 310-P2

- Measurement range:
 -50 °C ... +180 °C (-58 °F ... +356 °F)
- Needle: L = 130 mm, $\emptyset = 3 \text{ mm}$, blunt
- Cable: PTFE, L = 3 m

| Temperature | | Accuracy |
|-------------|-------------|----------|
| -50+60 °C | -58+140 °F | 0.6 °C |
| +60+180 °C | +140+356 °F | 0.9 °C |



TPX 310-P3

- Measurement range: -50 °C ... +180 °C (-58 °F ... +356 °F)
- Needle: L = 130 mm, Ø = 3 mm, blunt
- Cable: PTFE, L = 1 m

| Temperature | | Accuracy |
|-------------|-------------|----------|
| -50+60 °C | -58+140 °F | 0.5 °C |
| +60+180 °C | +140+356 °F | 0.8 °C |



TPX 310-P4

- · Measurement range: +100 °C ... +400 °C (+212 °F ... +752 °F)
- Needle: L = 50 mm, $\emptyset = 1.5 \text{ mm}$, blunt
- Cable: metal wrapped, L = 3 m, not waterproof

| Temperature | Accuracy |
|-------------------------|----------|
| +100+250 °C +212+482 °F | 1.1 °C |
| +250+400 °C +482+752 °F | 1.4 °C |



TPX 310-P5

- Measurement range: -50 °C ... +180 °C (-58 °F ... +356 °F)
 • Probe: L = 130 mm, Ø = 3 mm, blunt
- Cable: PTFE, L = 5 m

| Temperature | | Accuracy |
|-------------|-------------|----------|
| -5020 °C | -584 °F | 0.5 °C |
| -20+60 °C | -4+140 °F | 0.6 °C |
| +60+180 °C | +140+356 °F | 0.8 °C |
| | | |



TPX 310-P6

- Measurement range: -50 °C ... +180 °C (-58 °F ... +356 °F)
- Probe: L = 130 mm, Ø = 3 mm, blunt
- Cable: PTFE, L = 7.5 m

| Temperature | | Accuracy |
|-------------|-------------|----------|
| -50+60 °C | -58+140 °F | 0.7 °C |
| +60+180 °C | +140+356 °F | 1.0 °C |



TPX 310-P7

- Measurement range: -50 °C ... +180 °C (-58 °F ... +356 °F)
- Probe: L = 130 mm, Ø = 3 mm, blunt
- Cable: PTFE, L = 10 m

| Temperature | | Accuracy |
|-------------|-------------|----------|
| -50+60 °C | -58+140 °F | 0.9 °C |
| +60+180 °C | +140+356 °F | 1.1 °C |



How to order:

| Туре | Description | Part No. |
|-------------|---|------------|
| EBI 330-T30 | Single-use USB data logger, package unit: 10 pieces | 1340- 6332 |
| EBI 300 | Standard USB Data Logger | 1340-6330 |
| EBI 300 TE | USB Data logger with external temperature probe | 1340-6335 |
| TPC 300 | Replacement probe for EBI 300 TE | 1341-6331 |
| EBI 300 TH | USB Data logger with external humidity probe | 1340-6334 |
| TPH 400 | Replacement probe for EBI 300 TH | 1341-6336 |
| EBI 310 | High Precision USB Data Logger | 1340-6331 |
| EBI 310 TE | USB Data logger with external precision temperature probe | 1340-6337 |
| TPX 220 | Replacement probe for EBI 310 TE | 1341-6332 |
| EBI 310 TH | USB Data logger with external humidity probe | 1340-6336 |
| TPH 500 | Replacement probe for EBI 310 TH | 1341-6337 |
| EBI 310 TX | USB Data logger with temperature-two-channel-adapter | 1340-6339 |
| TPX 310 | Replacement adapter for EBI 310 TX | 1341-6335 |
| TPX 310-P1 | External sensor for EBI 310 TX | 1341-6338 |
| TPX 310-P2 | External sensor for EBI 310 TX | 1341-6339 |
| TPX 310-P3 | External sensor for EBI 310 TX | 1341-6340 |
| TPX 310-P4 | External sensor for EBI 310 TX | 1341-6341 |
| TPX 310-P5 | External sensor for EBI 310 TX | 1341-6342 |
| TPX 310-P6 | External sensor for EBI 310 TX | 1341-6343 |
| TPX 310-P7 | External sensor for EBI 310 TX | 1341-6344 |
| EBI 300-WM2 | Wall Mount for EBI 300 / 310 | 1340- 6341 |
| EBI 300 WM3 | Transportation mount for EBI 300/310 | 1340-6344 |









農園電 機台開發製作,電路和自動控制配電盤之設 FDA/USP/EPA美國藥典完全無菌等級減 器一處菌確放&冷鏈含儲溫歷度婚 對及製作,HMI/PLC 程式設計,系統規 菌/演畫是品異兩醇/變性乙醇。過乙酸,侧盤控。精密型/工業型電子天 劃, SCADA圖控監控軟體規劃整合設計。 過氣化氫。觀光型,例指釋過氧化氫溶 液(酚醛減菌劑/四級胺減菌劑/異丙醇 溶(兩) FLC可程式控制,SERVO 伺服馬 凌(酚醛減菌剂/四級胺減菌劑/異丙醇 強(兩) KERV 伺服馬 達、MS電磁閱關、荷重元、SENSOR酸鹼度/ 導電度/溶氧量/液(氣)體流量/壓力計。

美國原裝進口專業製藥工業ISO Class5 德國無線(有線)溫濕度壓力紀錄

11471 台北市内湖區新明路273巷6號1樓

1F., No.6, Ln. 273, Xinming Rd., Neihu Dist., Taipei City 11471, Taiwan (R.O.C.)

服務專線Tel: (02)8792-3722 服務傳真Fax: (02)8792-3761

電子信箱Email: info@grandever-biotech.com.tw 公司網址Website: www.grandever-biotech.com.tw





Xylem Analytics Germany Sales GmbH & Co. KG, ebro

Peringerstr. 10

85055 Ingolstadt, Germany

Tel +49 841 95478-0

Fax +49 841 95478-80

ebro@xyleminc.com

www.xylemanalytics.com