



Why you should be using a Refrigerator Thermometer

為什麼你應該使用冰箱溫度計

ebro TMX Series系列 - Digital Refrigerator Thermometer - 數位冰箱溫度計 Why do you need a refrigerator thermometer? The answer is quite simple - that is for you to

have complete visibility on the exact temperature inside refrigerators and freezers, which is especially critical for pharmaceutical and healthcare products.

為什麼需要冰箱溫度計?答案相當簡單 - 則讓您完全了解確切的冰箱和冰櫃內的溫度,這尤其重要用於製藥和保健產品。

All refrigerators used for pharmaceuticals, chemicals, greenhouses, blood banks, food and beverage must be monitored. Depending on the PIC/S definition, a deep freeze temperature should be kept below 15 °C and refrigerators should be within +2 °C to +8 °C. 所有用於製藥、化工、溫室血庫、食品和飲料的冰箱必須監測。取決於PIC/S 定義,深度冷凍溫度應保持在15 °C 以下而冰箱應在 +2 °C 至 +8 °C 內。



ebro TMX Series系列 - Digital Refrigerator Thermometer

- 數位冰箱溫度計

"...it is also important to know that the temperatures shown on the refrigerator display are not always representative of the entire fridge compartment."

...同樣重要的是要知道冰箱顯示屏上顯示的溫度並不總是代表整個冷藏室。"





Term 項目	Definition 定義
Temperature 溫度	Deep freeze 深度凍結: Below低於 -15 ℃
	In a refrigerator 在冰箱裡: +2 ℃ to +8 ℃
	Cold or Cool寒冷或冷凉: +8 ℃ to +15 ℃
	Room Temperature室溫: +15 °C to + 25 °C
	Ambient周圍的: The required storage temperature of non refrigerated medicinal product; usually stated on the product as "store below 25 °C" or "store below 30 °C". 非冷藏藥品所需的儲存溫度:通常在產品上註 明"儲存在 25 °C 以下"或"儲存在 30 °C 以下"。

Pic/s quide to good distribution practice for medicinal products Pic/s 藥品良好配送規範指南, pe 011-11 june 2014 https://picscheme.org/docview/3450

It is significantly important to use an appropriate refrigerator thermometer for fridges which are part of medical and healthcare processes that require precise temperature settings. Furthermore, we highly recommend to choose a thermometer with MAX/MIN display on it. When monitoring temperature sensible products, the minimum and maximum values must be recorded daily and documented manually. Last but not least, make sure to select a product with a calibration certificate. 對於需要精確溫 度設置的醫療和保健流程的一部分的冰箱,使用合適的冰箱溫度計非常重要。此外,我們強烈建議選擇 帶有 MAX/MIN 顯示的溫度計。在監測溫度敏感產品時,必須每天記錄最小值和最大值並手動記錄。最 後但並非最不重要的一點是,確保選擇帶有校准證書的產品。

Please take a look at ebro's TMX Digital Refrigerator Thermometer which offers some important features for you to consider when choosing your next refrigerator thermometer. 請査看 ebro 的 TMX 數字冰箱溫度計,它提供了一些重要功能,供您在選擇下一個冰箱溫度計時考慮。

ebro TMX Series系列- Digital Refrigerator Thermometer - 數位冰箱溫度計

- One channel or Two channels available—個管道或兩個管道可用
- Bottle insulates sensor option, prevent from rapid temperature changes when refrigerator door is opened瓶子絕緣傳感器選項,防止冰箱門打開時溫度急劇變化
- Wide temperature range -50 °C ~ +70 °C 寬溫度範圍 -50 °C ~ +70 °C
- Big display of current temperature and MAX/MIN value當前溫度和最大值/最小值 的大顯示屏
- 3 year calibration certificate 3 年校準證書
- Settable high and low limit 可設置的上下限
- Double alarm: Acoustic and Optical 雙重報警: 聲光報警
- Magnetic and pull-out stand attachment 磁性和拉出式支架附件



ebro TMX Series系列 - Digital Refrigerator Thermometer

- 數位冰箱溫度計



ebro TMX Series系列 - Digital Refrigerator Thermometer

- 數位冰箱溫度計





11471 台北市内湖區新明路273巷6號1樓

服務專線Tel: (02)8792-3722 服務傳真Fax: (02)8792-3761

電子信箱Email: info@grandever-biotech.com.tw 公司網址Website: www.grandever-biotech.com.tw

info.apac@xyleminc.com

www.xylem.com