# Annex 1 Practical Part- Webinar

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1. Screen / Presentation 屏幕/演示

2. Screen / Speaker 屏幕/揚聲器

The screenshot displays a Zoom webinar interface. On the left, a presentation slide is shown with the text: **-ebro-** – Second Salesmeeting, a xylem brand. Below this, a large black box contains the number **5**. The **-ebro-** logo is also visible in the bottom right corner of the slide. On the right, a blue screen displays the word **Speaker** in white text. A red double-headed arrow connects the two screens. At the bottom, a Zoom control bar is visible, featuring icons for **Chat**, **Hand heben** (Raise Hand), and **F&A** (Feedback). Two white arrows point from the **技術問題** (Technical Questions) label on the left and the **問題與解答** (Questions and Answers) label on the right towards the **Chat** and **F&A** icons respectively. A red **Verlassen** (Leave) button is located in the bottom right corner of the control bar.

# Information – LIVE Webinar via Zoom

1. Screen / Presentation

2. Screen / Speaker

The screenshot displays a Zoom webinar interface. On the left, a screen share is active, showing a blue background with the **-ebro-** logo (a xylem brand) and the text "Second Salesmeeting". A large white number "5" is centered on a black rectangular area. A red double-headed arrow points from this screen to the right, where a "Speaker" view is shown against a blue background. Below the main content, a text overlay reads "Speaker view can be swapped with a double click". At the bottom, the Zoom control bar is visible, showing the name "Robert Streller", "Audio-Einstellungen", "Chat", "Hand heben", "F&A", and a "Verlassen" button.

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# Brief Introduction

# Webinar Speaker:

## Sebastian Schwarz - Market Manager Pharma 市場經理製藥

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- 8 years experience as a consultant in a GMP controlled environment 8 年在 GMP 受控環境中擔任顧問的經驗
- Projects in the areas of calibration, qualification and validation for pharmaceutical manufacturers throughout Europe 為整個歐洲的製藥商提供校準、鑑定和驗證領域的項目

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# - ebro - Markets, Applications And Devices 市場、應用和設備

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# Regulatory 監管

**Regulatory bodies are organisations set up by 監管機構**是由以下機構設立的組織

Government with responsibility to monitor, guide and control various industry sectors in the interests of protecting consumers – *Protect, Promote and Maintain the health and safety of the public.*

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**Widely Known Regulatory Authorities 廣為人知的監管機構**

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EUROPEAN MEDICINES AGENCY  
SCIENCE · MEDICINES · HEALTH

Defence and  
Security  
Accelerator

***\*There are other Organizations and Associations in the Industry that influence the Pharmaceutical Industry very highly! (ex. The PDA, ISPE, ABPI, etc.) 工業內還有其他組織和協會對製藥工業影響很大！ (例如 PDA、ISPE、ABPI 等)***

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**However, there is a different approach:** 但是，有一種不同的方法：

- Annex 1 附件 1：  
Biological indicators should be considered as an **additional method** for monitoring the sterilization. 應考慮將生物指示劑作為監測滅菌的附加方法。
- FDA美國食品和藥物管理局：  
In general, the **biological indicator should be placed** adjacent to the temperature sensor to assess the correlation between microbial lethality and predicted lethality based on thermal input. 一般來說，生物指示劑應放置在溫度傳感器附近，以評估微生物致死率和基於熱輸入預測的致死率之間的相關性。

# Important norms regarding sterilization processes

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What kind of sterilization processes have you used before? 您以前使用過什麼樣的滅菌流程？

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# 4


# Validation of sterilization processes

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# Sterilization processes 滅菌過程


## Healthcare / Pharmaceutical 醫療保健/製藥

1. Moist Heat Sterilization 濕熱滅菌
2. Dry Heat Sterilization / Depyrogenation 乾熱滅菌/除熱原
3. Ovens / Vacuum Ovens 烘箱/真空烘箱
4. H<sub>2</sub>O<sub>2</sub> Sterilization 過氧化氫滅菌
5. EtO Sterilization 環氧乙烷滅菌

# The V-Model of qualification according to EU GMP Annex 15

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# Conceptual Design – User Requirement Specification

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### Knowing your risks, minimize them during conceptual design!

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**Check of programs (with empty chamber):** 檢查程序（空腔）：

- Fractionated pre-vacuum 分段預真空
- Heating phase 加熱階段
- Equilibration time 平衡時間
- Sterilization time 滅菌時間
- Drying and Cooling 乾燥和冷卻

# Qualification Phase - PQ 資格階段 - PQ

**PQ** – Performance qualification of all procedures & sequences: 所有程序和序列的性能鑑定：

Checks 檢查:


**Conductivity of the feed water** 給水的電導度

CT 830 Set Conductivity Meter 電導度套組

# Qualification Phase 資格階段 - PQ

## Vacuum test 真空測試


## Ebro Solutions 解決方案

- EBI 12 / EBI 11
- Temperature – Pressure data logger 溫度壓力記錄器

# Qualification Phase 資格階段 - PQ

## Bowie Dick-test 鮑伊迪克測試


Electronic Bowie Dick Test according to DIN EN ISO 11140-4 符合 DIN EN ISO 11140-4 的電子鮑伊迪克測試

EBI 16 Bowie Dick Test Data Logger 鮑伊迪克測試記錄器

# Qualification Phase 資格階段- PQ

## Empty chamber profile 空腔輪廓

Goal is to map the autoclave to find possible cold spots  
目標是測繪高壓釜以找到可能的冷點

### Ebro Solutions 解決方案

- EBI 12 / EBI 11
- Temperature -data logger 溫度記錄器

# Moist Heat Sterilization Validation 濕熱滅菌驗證

## 3 successful validation runs Thermometric check with full loading 3次成功的運行滿載的溫度檢查驗證

- One probe attached to the Primary Reference 一個探頭連接到主要參考
- One probe attached or near the RTD - Controller 一個探頭連接或靠近標準件 - 控制器
- Two probes in containers filled with the testing medium 兩個探針置於裝有測試介質的容器中的
- A minimum of 3 probes each located in separate areas of each crate or tray 至少 3 個探頭，每個探頭位於每個板條箱或託盤的不同區域
- Initial Temperature and Data Collection Points 初始溫度和數據採集點
- A schematic drawing to show the placement of all probes 顯示所有探頭位置的示意圖
- Defined loading patterns 定義的裝載模式

### Ebro Solutions 解決方案

- EBI 12 / EBI 11
- Temperature – Pressure data logger 溫度-壓力記錄器

# Reaching the sterilization conditions 達到滅菌條件

Recommended minimum number of temperature test locations (sensors)

Base: about 12 Measuring points/1000 l usable volume.

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# Reaching the sterilization conditions 達到滅菌條件

The existence of saturated steam in the usable area and within the loading can be considered as reached, **if all temperatures measured** in the usable area and within the loading during the hold time: 如果在保持時間內在可用區域和裝載內測量到的所有溫度，則可以認為在可用區域和裝載內存在飽和蒸汽：

- are not below the sterilization temperature 不低於滅菌溫度
- are not more than 3K (large sterilizer and small since 2016) 不超過 3K (自 2016 年以來大型滅菌器和小型)
- do not diverge by more than 2K 偏差不要超過 2K
- equilibration time 15 s to 800 liter, 30 s for larger sterilizers 平衡時間 15 秒至 800 公升，大型滅菌器為 30 秒
- minimum hold times 121 °C for 15 min; 126 °C for 10 min; 134 °C for 3 min 最短保持時間 121 °C 保持 15 分鐘； 126°C 10 分鐘； 134°C 3 分鐘
- $F_0$ -Value minimum 15 min  $F_0$  值至少 15 分鐘

The saturated steam temperature that is calculated with the measured pressure is to be considered as measured temperature. 以測量壓力計算的飽和蒸汽溫度作為測量溫度。

# Measuring values outside the tolerances 測量值超出公差

What measures need to be taken in case the results do not correspond to the specifications? 如果結果不符合規範，需要採取哪些措施？

1. Check if the loading is correct? 檢查裝載是否正確？
2. Positioning of the measuring sensors 測量傳感器的定位
3. Repeated calibration of the measuring system 測量系統的重複校準
4. Adjust process parameters by manufacturer or service company (e.g. confining / chamber pressure, pressure period) 由製造商或服務公司調整流程參數（如圍壓/腔壓、壓力週期）
5. Alternative medical products? (materials, construction) 替代醫療產品？（材料、結構）

# Depyrogenation 去熱原

## Batch Hot Air Ovens 批次熱風烘箱

- The glass containers are loaded into batch ovens and then sterilized generally 2-6 hours depending on the temperature. 將玻璃容器裝入批次烘箱，然後根據溫度進行消毒，一般需要 2-6 小時。
- Batch ovens typically require 12 to 24 measuring points depending on size. 批次烘箱通常需要 12 到 24 個測量點，具體取決於尺寸。
- Use of the EBI TIB Thermal Insulation Box is necessary when using either the EBI 12-T22X or the EBI 12-T421. 使用 EBI 12-T22X 或 EBI 12-T421 時，必須使用 EBI TIB 隔熱箱。
- The EBI 40 is a practical alternative for the Batch Ovens EBI 40 是批次烘箱的實用替代品

# Depyrogenation 去熱原

## Continuous Hot Air Tunnels 連續熱風隧道

- The glass containers are run via conveyors through the tunnel, typically at 300°C or higher, from 15 to 30 minutes. 玻璃容器通過傳送帶通過隧道，通常在 300°C 或更高溫度下運行 15 到 30 分鐘。
- Tunnels usually require 3 rows of 5 monitoring points about 3 feet apart, totally 15 measuring points 隧道通常需要3排5個相距約3英尺的監測點，共15個測量點
- EBI 12-T421 is recommended with the EBI TIB Thermal Insulation Box. 建議將 EBI 12-T421 與 EBI TIB 隔熱箱一起使用。

## Ovens / Vacuum Ovens 烘箱/真空烘箱

- Vacuum Ovens are used for drying heat sensitive materials such as powder to extract moisture 真空烘箱用於乾燥熱敏材料例如粉末以粹取水分
- The vacuum oven enables drying treatment at a lower temperature for specimens that cannot be treated by conventional high-temperature drying. • 真空烘箱可以在較低溫度下對常規高溫乾燥無法處理的樣品進行乾燥處理。
- Validation is necessary to avoid too fast, too slow or unevenly drying of the goods 必須進行驗證以避免貨物太快、太慢或乾燥不均勻
- Typically, temperature testing, evaporation, drying 通常，溫度測試、蒸發、乾燥

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# Survey

Which information and reports are useful for your process

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**-ebro-**  
a xylem brand

# 7

## Evaluation of a steam sterilization process with the Winlog.validation 使用 Winlog.validation 評估蒸汽滅菌過程

# Requirements for the measuring technology as per EN 285 / ISO 17665 根據 EN 285 / ISO 17665 對測量技術的要求

## Sensors for temperature measurement 溫度測量傳感器

- Sensor diameter 傳感器直徑 <  $\text{Ø } 2 \text{ mm}$  ( $3,1\text{mm}^2$ )
- Sensor length 傳感器長度 given by application
- Range 範圍  $0 \dots 140^\circ\text{C}$
- Accuracy 精準度  $\pm 0,25\text{K}$

## Sensors for pressure measurement 壓力測量傳感器

- Range 範圍 <  $50 \dots 4000\text{mbar}$   
( $< 5\text{kPa} \dots 400\text{kPa}$ )
- Accuracy 精準度  $\pm 20 \text{ mbar at } 25^\circ\text{C}$
- Connector 連接器 ( $\pm 2 \text{ kPa at } 25^\circ\text{C}$ ) no connector required (chamber pressure) 無需連接器 (腔室壓力)

# Validation of the thermal process in sterilizers 驗證滅菌器中的熱過程

## Suitable temperature logger 合適的溫度記錄器

EBI 12 temperature logger with 1.2 to 1.5 mm flex probes radial and axial

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EBI 11 temperature logger with 2 mm rigid probes (diameter 1.95 mm)

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## Suitable pressure logger 合適的壓力記錄器

EBI 12 pressure logger with 1.5 mm flex probe with luerlock

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EBI 11 pressure logger with and without luerlock

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Validated by independent company like TÜV  
經 TÜV 等獨立公司驗證


Calibration certificates 校準證書

IQ, OQ documentation 確效文件


# Validation evaluation 驗證評估

What would a typical steam sterilization report look like? 典型的蒸汽滅菌報告是什麼樣的？

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
$$F_0 = \Delta t \sum 10^{\frac{T-121.1}{z}}$$


Δt – the time interval between two temperature readings 兩次溫度讀數之間的時間間隔

T – the temperature at time t of the product under sterilization 滅菌產品在時間 t 的溫度


z – temperature coefficient (assumed as 10 C) 溫度係數（假設為 10°C）

# Validation evaluation 驗證評估

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# Validation evaluation 驗證評估

Result

**Passed 通過**

Criterion	Nominal value	Actual value
✓ Sterilization temperature band	134.00 - 137.00 °C	134.06 - 135.55 °C
✓ Maximum equilibration time	≤ 15 Sec.	1 Sec.
✓ Minimum sterilization time	≥ 180 Sec.	337 Sec.
✓ Lethality	≥ 15.00	148.86
✓ #15114714 [1]	≥ 15.00	150.85
✓ #15114718 [1]	≥ 15.00	148.86
✓ #15114712 [1]	≥ 15.00	150.98
✓ #15114744 [2]	≥ 15.00	154.16
✓ #15114671 [1]	≥ 15.00	151.59
✓ #15114717 [1]	≥ 15.00	149.83
✓ #15114715 [1]	≥ 15.00	152.18
✓ Variance	≤ 2.00 K	0.31 K
ⓘ Pressure velocity decrease (evacu...)		133.8 kPa/Min.
ⓘ Pressure velocity decrease (drying)		293.5 kPa/Min.
ⓘ Pressure velocity increase (equilibr...)		102.2 kPa/Min.

Ok    Set manually to 'passed'

Further questions? Don't hesitate to contact me under: 進一步的問題？  
請隨時通過以下方式與我聯繫：

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## Webinar

### Stale air is so last year!

With the CO<sub>2</sub>-measuring device RM 100 for air quality measurement

Tuesday, **29th March**

**2:00 p.m. to 2:30 p.m.**

Roland Trübswetter

Productmanager Food

**book your seat!**

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GRANDEVER BIOTECHNOLOGY CO., LTD.