

Professional Data Loggers for Medical Industry 應用於醫療行業的專業數據記錄器

Validation and Routine monitoring of steam sterilizer and washer-disinfector processes in hospitals 蒸汽滅菌器的驗證和日常監測
醫院的清洗消毒流程



Agenda

議程

1. Company Introduction 公司介紹
2. Where and when do you need to use data logger for medical industry 醫療行業何時何地需要使用數據記錄儀
3. Validation Loggers for Steam Sterilizers Process 蒸汽滅菌器過程的驗證記錄器
4. Validation Loggers for Washer-Disinfectors Process 清洗消毒器過程的驗證記錄器
5. Success Stories 成功案例
6. Solutions for Other Devices 其他設備的解決方案
7. Q&A 問答

1. Company Introduction

公司介紹

Introduction 介紹



Iven Kruse

Introduction 介紹



Helmut Schmidt

Ebro Today

- Founded in Freiburg Germany in 1968
1968 年在德國弗萊堡成立
- Focused on Medical, Pharma and Lab/
Industrial markets 專注於醫療、製藥和實驗室/工業市場
- Global coverage with over 100 distributors
worldwide 全球覆蓋超過 100 家分銷商
- Xylem Analytics Germany Sales GmbH &
Co. KG (XA Sales) since July 2016
自 2016 年 7 月起加入 Xylem Analytics
Germany Sales GmbH & Co. KG (XA Sales)



2. Who? Where? When? Why? 誰？在哪裡？ 什麼時候？為什麼？

WHO needs the Professional Data Logger ? 誰需要專業的數據記錄器?



WHERE? 在哪裡?

Sterile Goods Cycle

無菌物品循環



Operating Room Doctor 手術室醫生



WHERE can you use our data logger? 在哪裡您可以

使用我們的數據記錄器?

Process Validation & Routine monitoring for 流程驗證和例行監控

Washer-Disinfector 清洗消毒器

- Washer disinfector process 清洗消毒器消毒程序
- Washer disinfector for endoscopes processes 用於內窺鏡過程的清洗消毒器
- Human waste container washer disinfector 人體廢物容器清洗消毒器



Steam Sterilizer 蒸汽滅菌器

- Electronic B&D test 電子式Bowie Dick蒸汽滲透真空測試
- Steam sterilizer processes 蒸汽滅菌器過程



Others 其他

- DAC processes 牙科高壓滅菌器/ Careclave 牙科器械的卡盤護理 (Dental牙科)
- Plasma sterilizer processes (H₂O₂) 電漿態過氧化氫滅菌流程
- Ethylene Oxide processes (EtO) 環氧乙烷流程
- Formaldehyde processes (HCHO) 甲醛流程



WHY? 為什麼?

U.S.A 美國

- Biological Monitoring 生物的監測
- Chemical Monitoring 化學的監測

Guidelines: 準則

ANSI/AAMI ST79

European, Australia 歐洲, 澳洲


- Validation 驗證
- Physical Monitoring 物理的監測

Guidelines: 準則

EN 285
ISO 17665

3. Validation for Steam Sterilizer Process 蒸汽滅菌 器流程驗證

Reprocessing of medical devices 醫療器械再處理



High standards are required
需要高標準:

- Protection of patients 保護患者
- Protection of employees 保護職員
- Legal security 法規保障
- Savings capacity 儲積能力

Definition „Validation“ 定義“驗證”

What is a validation 什麼是驗證？

Repeatable verification of the fact, that a process permanently reaches the required results. 事實的重複驗證是一個永久達到需求結果的過程。

Where or when does a validation have to be performed 必須在何時何地進行驗證？

In any case where the positive result of a process cannot be verified by simply one measurement. 在任何情況下，僅通過一次測量無法驗證過程的正面結果。

Definition „Validation“ 定義“驗證”

- A validation includes the Installation Qualification (IQ), the Operation Qualification (OQ) and the Performance Qualification (PQ). 驗證包括安裝認證 (IQ)、操作認證 (OQ) 和性能認證 (PQ)。
- A validation is clear prove for the fact, that procedures, processes, equipment, materials, process steps or systems actually lead to the expected results. 驗證清楚地證明程序、過程、設備、材料、流程步驟或系統實際上引導預期結果的事實。
- The results are summarized and annotated in a validation report, which helps to work out and realize propositions for process optimization. 結果在驗證報告中進行了總結和註釋，這有助於制定和實現流程優化的建議。



Legal basis and standards 法律依據和標準

Standards and guidelines for steam sterilizer processes: 蒸汽滅菌器流程的標準和指南：

DIN EN 285 large sterilizers 大型滅菌器

EN 13060 small steam sterilizers 小型蒸汽滅菌器 (< 60 liters公升)

EN ISO 17665 sterilization of products for health care – humid heat -requirements for the development, validation and control of the application of a sterilization process for medical products
醫療保健產品的滅菌 - 濕熱 - 醫療產品滅菌過程應用的開發、驗證和控制要求


EN ISO 11140-4 Routine control of steam sterilizers - alternative Bowie Dick Tests 蒸汽滅菌器的常規控制 - 電子式Bowie Dick蒸汽滲透真空替代測試

Typical annual re-qualification is necessary 典型的年度重新認證是必要的(ISO 17665-2 part 12.4).

Installation Qualification 安裝認證 (IQ)

1. Summary of the results 結果總結
2. Deviations and defects 偏差和缺失
3. General information 一般信息
4. Task formulation 任務制定
5. Identification of the device 設備識別
6. Check the documentation for completeness 檢查文檔的完整性
7. Organisational preconditions for the validation 驗證的組織前提條件
8. Installation control: environmental conditions, operating material, equipment, maintenance 安裝控制：環境條件、操作材料、設備、維護
9. Description of the general state of the sterilizer 滅菌器一般狀態說明
10. Company organisation – control of the quality assurance measures (trainings, daily operation, work instructions) 公司組織-品質保證措施的控制（培訓、日常操作、工作指導）
11. Function 功能
12. Used measuring equipment 使用的測量設備

Prüfpunkte	Entsprechen: Ja / Nein	Bemerkung / Kommentar	Datum / Unterschrift
Allgemeine Angaben			
Allgemeine Angaben vollständig			
Aufgabenstellung vollständig			
Identifikation des Gerätes vollständig			
Prüfung der Dokumentation vollständig			
Organisatorische Vorbedingungen zur Validierung			
Installationskontrolle			
Umgebungsbedingungen			
Betriebsmittel			
Ausrüstung			
Allgemein			
Ausführung und Konstruktion			
Wartung			
Beschreibung des Allgemeinzustandes des Sterilisators			
Betriebsorganisation			
Arbeitsanweisungen			
Qualitätssicherung des täglichen Betriebs			
Schulungen			
Funktion			
Verwendete Messmittel			



Operation Qualification 操作認證(OQ)

1. General information 一般信息
2. Summary of the results and deviations / defects 結果和偏差/缺失的總結
3. Preconditions for the Operation Qualification 操作認證的前置條件
 - Collect medical products for the reference loading 收集醫療產品以供參考裝載
 - Record medical products in packing lists 在裝箱單中記錄醫療產品
 - Record the test configuration 記錄測試配置
4. Checks 檢查
 - Conductivity of the feed water 給水的電導率
 - Vacuum test 真空測試
 - BD-test and/or Helix-test 抽真空BD 測試和/或 Helix 測試
 - Empty chamber profile 空腔輪廓
 - Thermometric check with partial loading 部分加載的測溫檢查
 - Thermometric check with full loading 滿載溫度測量檢查



Performance Qualification 性能認證(PQ)

1. General information 一般信息
2. Summary of the results and deviations / defects 結果和偏差/缺失的總結
3. Preconditions for the Performance Qualification 性能認證的前提條件
4. Checks 檢查
 - BD-test 蒸汽滲透真空測試
 - **Three** times product check 三次產品檢查
 - Proof of product compatibility 產品兼容性證明
 - Proof of protection from recontamination 防止再污染的證明
 - Proof of the identification ability 識別能力證明
5. Determination of routine controls 常規控制的量測



Validation of processes in steam sterilizers

蒸汽滅菌器中的驗證流程



Thermometric check of the **partial loading** using dataloggers
使用數據記錄器對**部分加載**進行溫度測量

Validation of processes in steam sterilizers

蒸汽滅菌器中的驗證流程




Thermometric check of the **full loading** using dataloggers
使用數據記錄器對**滿載**進行溫度測量



Interpretation of the thermometric measurement result

溫度測量結果的解讀



Reaching the sterilization conditions 達到滅菌條件

The existence of **saturated steam** in the usable area and within the loading can be considered as reached, if all temperatures measured in the usable area and within the loading during the hold time: 如果在保持時間內在可用區域和裝載內測量到的所有溫度，則可以認為在可用區域和裝載內存在**飽和蒸汽**：

- are **not below** the sterilization temperature **不低於**滅菌溫度
- are not more than **3K** above the sterilization temperature 高於滅菌溫度不超過 **3K**
- do not diverge by more than **2° C** 偏差不要超過 **2° C**
- equilibration time **15 s** to 800 liter, **30 s** for larger sterilizers 平衡時間 **15 秒**至 800 升，大型滅菌器為 **30 秒**
- **minimum** hold times 121° C for 15 min; 126° C for 10 min; 132° C for 3,5 min 134° C for 3 min
最短保持時間 121° C 保持 15 分鐘； 126° C 10 分鐘； 132° C 3.5 分鐘 134° C 3 分鐘
- F₀-Value minimum **15 min** F0 值至少 **15 分鐘**

The **saturated steam** temperature that is calculated with the measured pressure is to be considered as measured temperature. 以實測壓力計算的**飽和蒸汽**溫度作為實測溫度。



What if measuring values outside the tolerances 如果測量值超出公差範圍怎麼辦

What measures need to be taken in case the results **DO NOT** correspond to the specifications? 如果結果不符合規格，需要採取哪些措施？

1. Check if the loading is correct? 檢查加載是否正確?
2. Positioning of the measuring sensors 測量傳感器的定位
3. Repeated calibration of the measuring system 測量系統的反復校準
4. Adjust process parameters by manufacturer or service company (e.g. confining / chamber pressure, pressure period)
由製造商或服務公司調整流程參數（如限制/腔壓、壓力週期）
5. Alternative medical products? (materials, construction) 替代醫療產品？（材料、結構）



Requirements for the measuring technology as per EN 285 / ISO 17665

根據 EN 285 / ISO 17665 對測量技術的要求

Temperature data logger 溫度數據記錄器

- Temperature sensor made from platinum resistors **Pt class A** as per IEC 60751 or thermocouples class 1 as per IEC 60584-2, 由符合 IEC 60751 的 **A 級鉑電阻器** 或符合 IEC 60584-2 的 1 級熱電偶製成的溫度傳感器
- Diameter of the sensor \varnothing less **than 2 mm** 傳感器直徑小於 **2 mm**
- Temperature measuring range **0 to 150 °C** 溫度測量範圍 **0 至 150 °C**
- Measuring interval 測量間隔 ≤ 2.5 s
- Resolution 解析度 0.1 K
- Accuracy of the temperature sensor at least 溫度傳感器的精度至少 0.5 K
- Accuracy of the temperature recording device 溫度記錄裝置的精度 ± 0.25 K
- Calibration of the temperature measuring technology according to manufacturer instructions; a temperature that lies within the sterilization temperature must be calibrated 根據製造商說明校準溫度測量技術；必須校準處於滅菌溫度範圍內的溫度 (DIN EN 285 Pt. 26.4.9)



Requirements for the measuring technology as per EN 285 / ISO 17665

根據 EN 285 / ISO 17665 對測量技術的要求

Pressure Data Logger 壓力數據記錄器:

- Temperature-compensated pressure sensor 溫度補償壓力傳感器
- Measuring rate 測量率 ≤ 1 s
- Measuring range of the pressure measuring device 壓力測量裝置的測量範圍: **0 to 4 bar** (0 to 400 kPa)
- Resolution 解析度: 0.04 bar (4 kPa)
- Accuracy of the pressure measuring device $\pm 0.5\%$ of the measuring range, as long as the environment temperature is $(20 \pm 3^\circ\text{C})$ 壓力測量裝置的精度為測量範圍的 $\pm 0.5\%$ ，只要環境溫度為 $(20 \pm 3^\circ\text{C})$
- The device must have a valid certificate 設備必須具有有效證書
- Calibration according to manufacturer instructions, traceable to a national standard (e.g. DKD), (DIN EN 285 Pt. 26.5) 根據製造商說明進行校準，可溯源至國家標準（例如 DKD），（DIN EN 285 Pt. 26.5）



Validation report for sterilization processes

滅菌流程驗證報告

Validering Winlog.validation

Descr.: Tandhälsdan Tyresö 2017-11-30 Temp.prof.
 30.11.2017 10:55:59 Winlog.validation 3.71

Stockholm

General parameters

Device SOP
 Stern Weber 17+ -

Program Standard
 Temp.profil DIN EN 13060

Creator Charge number
 Admin -

Responsible Validated
 Administrator 08.03.2018 15:07:16

Remark
 s/n: 17GP0028

Loggers in process


# 15177422	Calibration date 28.09.2017 09:48:00
# 15187298	Calibration date 28.09.2017 09:46:04
# 15187311	Calibration date 28.09.2017 09:47:36
# 15194093	Calibration date 28.09.2017 09:46:56
# 15202146	Calibration date 28.09.2017 09:46:39
# 15232080	Calibration date 28.09.2017 09:47:10

Overall result

Passed 已通過

Detailed results	Nominal	Actual
Pressure velocity increase (equilibration)		76.7 kPa/Min.
Pressure velocity decrease (drying)		296.2 kPa/Min.
Pressure velocity decrease (evacuation)		146.7 kPa/Min.
Max. variance	<= 2.00 K	0.24 K
Min. sterilization time	>= 180 sec.	270 sec.
Max. equilibration time	<= 15 sec.	5 sec.
Temperature range	134.00 - 138.00 °C	134.81 - 135.18 °C

Total measurement



Executed by

Verified by

Kwalifikacja procesowa plastica

Descr.: 1st_16_11_2021-AC4-P1
 05.01.2022 13:53:17 Winlog.validation 3.71

Warsaw

General parameters

Device SOP
 AC4 -

Program Standard
 P1 -

Creator Charge number
 Admin 3

Responsible Validated
 Admin 13.04.2022 10:29:40

Loggers in process


# 19705712	Calibration date 02.03.2021 08:04:25
# 19706213	Calibration date 13.01.2021 06:50:03
# 19706214	Calibration date 13.01.2021 06:49:09
# 19706215	Calibration date 13.01.2021 06:48:05
# 19706216	Calibration date 13.01.2021 06:47:48
# 19706217	Calibration date 13.01.2021 06:50:14

Overall result

Passed 已通過

Detailed results	Nominal	Actual
Pressure velocity increase (equilibration)		23.1 kPa/Min.
Pressure velocity decrease (drying)		107.8 kPa/Min.
Pressure velocity decrease (evacuation)		192.7 kPa/Min.
Max. variance	<= 2.00 K	1.16 K
Min. sterilization time	>= 330 sec.	357 sec.
Max. equilibration time	<= 15 sec.	6 sec.
Temperature range	134.00 - 137.00 °C	134.01 - 135.87 °C

Total measurement




Executed by Administrator

Verified by

What is Bowie Dick Test? 什麼是Bowie Dick蒸汽滲透真空測試?

To check if the air removal and steam penetration are successful
檢查排氣和蒸汽滲透是否成功

Original Test Pack 原始測試包



Alternative BD Test 替代BD測試—
Electronic Bowie-Dick Test
電子式Bowie Dick蒸汽滲透真空測試

Alternative BD Test 替代BD測試—
Chemical Indicator 化學指示劑




Alternative BD Test 替代BD測試—
Chemical Indicator 化學指示劑



Definition „Routine Control“ Daily Bowie Dick Test 定義“常規控制”每日鮑伊迪克測試


- **Routine control** is the periodical verification to **check if the operating performance** of the washer-disinfector or sterilizer meets the limit values that were defined during validation. 例行控制是定期驗證，以檢查清洗消毒器或滅菌器的**操作性能**是否滿足驗證期間定義的限值。



Electronic Bowie & Dick Test 鮑伊迪克測試 EBI 16



- **Electronic BD** test according to EN 285 / ISO 17665
根據 EN 285 / ISO 17665 進行**電子 BD** 測試
- Air removal steam penetration test to verify the quality of sterilization
除氣蒸汽滲透測試，驗證滅菌質量



4. Validation for Washer-Disinfector Process 清洗消毒器 流程的驗證

Legal Background法規背景

- EN ISO 15883
routine monitoring and validation of washer disinfectors 清洗消毒器的常規監測和驗證
- EN ISO 15883-1:
General requirements, definitions and verifications 一般要求、定義和驗證
- EN ISO 15883-2:
Requirements and checks for washer-disinfector devices with thermal disinfection for surgical instruments, anaesthesia instruments etc. 對手術器械、麻醉器械等進行熱消毒的清洗消毒裝置的要求和檢查。
- EN ISO 15883-3:
Requirements and checks for washer-disinfector devices with thermal disinfection for containers for human excretions 人體排泄物容器熱消毒清洗消毒裝置的要求和檢查
- EN ISO 15883-4:
Requirements and checks for washer-disinfector devices with chemical disinfection for thermolabile endoscopes 對耐熱內窺鏡進行化學消毒的清洗消毒裝置的要求和檢查
- WFHSS Guideline No. 02: Check List for Procurement of Medical Devices pursuant to EN ISO 17664:2004
WFHSS 指南第 02 號：根據 EN ISO 17664:2004 採購醫療器械的清單

Washer disinfectors - Validation

清洗消毒器 - 驗證

Influences to the automatic cleaning and disinfection
對自動清洗消毒的影響

1. Washer disinfecter device 清洗消毒裝置:

- Cleaning programs 清潔程序
- Disinfection programs 消毒計劃
- Rinse water circulation 沖洗水循環

2. Loading 裝載:

- Arrangement towards the rinsing mechanics 對沖洗機制的安排
- Rinsing dead zones 沖洗死區
- Hollow spaces 空洞空間



Washer disinfectors - Validation

清洗消毒器 - 驗證

Influences to the automatic cleaning and disinfection 對自動清洗消毒的影響

3. Equipment設備:

- Water 水
- Steam 蒸汽
- Compressed air 壓縮的空氣
- Energy supply 能源供應

4. Process Chemicals: 流程化學品

- Detergents 洗滌劑
- Disinfectants 消毒劑
- Neutralization agents 中和劑
- Rinsing agents 沖洗劑



Washer disinfectors - Validation

清洗消毒器 - 驗證

Influences to the automatic cleaning and disinfection

對自動清洗消毒的影響

5. Instrument Design 儀器設計:

- Materials 材料
- Joints 連接處
- Hollow spaces 空洞空間

6. Load Carrier 裝載載體:

- Rotating Arms 旋轉臂
- Instrument holders 儀器支撐
- Spray nozzles 噴嘴




Washer disinfectors - Validation

清洗消毒器 - 驗證

What is essential for the validation 驗證的必要條件是什麼？

- Temperature increase and hold time during the cleaning step and thermal disinfection 清潔步驟和熱消毒過程中的溫度升高和保持時間
- Pressure changes 壓力變化
- pH value 酸鹼度值
- Conductivity effect of the demineralized water 軟化水的電導度效應
- Conductivity last rinsing water before draining 排水前最後沖洗水的電導度



Washer disinfectors - Validation

清洗消毒器 - 驗證

The validation consists of 驗證包括：

- Installation Qualification 安裝認證 IQ
- Operation Qualification 操作認證 OQ
- Performance Qualification 性能認證 PQ



Washer disinfectors - Validation

清洗消毒器 - 驗證

Installation Qualification 安裝認證 IQ

The Installation Qualification confirms that 安裝認證確認

- The WD and all ist accessories have been delivered properly and the installation has been performed professionally. WD 和所有ist 附件均已正確交付，並且安裝已專業完成。
- The connection to the matching supplies, water, power supply and process chemicals dosing pump system have also been completed perfectly. 與配套供應、水、電力供應和流程化學品計量泵系統的連接也已完美完成。



Washer disinfectors - Validation


清洗消毒器 - 驗證

Operation Qualification 操作認證OQ

The Operation Qualification confirms that 操作認證確認

The complete WD unit and all the supplies are **in accordance with the manufacturers requests and** match the demands laid down in the **DIN EN ISO 15883**.

完整的WD裝置和所有供應品**均符合製造商的要求**，並符合**DIN EN ISO 15883** 中規定的要求。



Washer disinfectors - Validation

清洗消毒器 - 驗證

Cleaning and Disinfection Performance 清潔消毒性能

Checking the cleaning performance is done in two steps which can be combined, if necessary: 檢查清潔性能分兩個步驟進行，如有必要，可以合併：


1. Test instruments (test body) with defined test soiling 具有定義測試污染的測試儀器（測試體）
2. Typical instruments, contaminated by use 典型儀器，因使用而污染



Washer disinfectors - Validation

清洗消毒器 - 驗證

Test instrument clamp according to Crile
with 0.1mg heparinised sheep blood with Protamine sulfate
根據 Crile 以0.1mg 肝素化羊血和硫酸魚精蛋白測試儀器夾具



Washer disinfectors - Validation

清洗消毒器 - 驗證

Test body and reference loading
with contaminated, general instruments
以使用受污染的通用儀器測試體和參考裝載



Washer disinfectors - Validation

清洗消毒器 - 驗證

Step 步驟 1 - Check optical cleanliness 檢查光學清潔度

The test instruments
(clamps according to Crile)
as well as the marked
contaminated instruments
are checked
regarding optical
cleanliness.

關於光學清潔度（根據 Crile
的夾具）測試儀器和標記受
污染儀器進行檢查



Washer disinfectors - Validation

清洗消毒器 - 驗證

Step 步驟 2 - Check for protein remains using the Biuret method 使用雙縮脲法檢查蛋白質殘留

- Remove the protein remains with SDS solution (sodium dodecyl sulfate; 1% in water) which is included in the test set or available in your own hospitals pharmacy 使用 SDS 溶液（十二烷基硫酸鈉；1% 的水溶液）去除蛋白質殘留物，該溶液包含在測試集中或您自己的醫院藥房有供應




Washer disinfectors - Validation

清洗消毒器 - 驗證

Biuret test set by Miele 德國美諾 Miele 雙縮脲測試儀



Biuret method with the swab procedure protect
雙縮脲法與擦拭程序保護



Washer disinfectors - Validation

清洗消毒器 - 驗證

Pressure measurement 壓力測量

Sufficient water pressure is precondition for the functioning of the spraying system 用於噴灑系統的功能前提是水壓充足

- A pressure loss already occurs due to the blood proteins rinsed off and dissolved from the medical products. 由於血液中的蛋白質被沖洗掉並從醫療產品溶解導致壓力損失。
- Pre-cleaning detergents can cause a major pressure drop. 預清潔劑可能會導致一個主要的壓力下降。
- Pressure changes affect the rotation of the spray arms. 壓力變化影響噴霧臂的旋轉。
- **By using pressure loggers which record pressure and temperature in relation to time permanently during the complete process, these influences can be measured and evaluated.** 通過使用在整個過程中永久記錄與時間相關的壓力和溫度的壓力記錄器，可以測量和評估這些影響。



Washer disinfectors - Validation

清洗消毒器 - 驗證

An accident set should contain the following:
事故套裝應包含以下內容：



Description 描述

- A simple possibility of batch control is offered by the following product compilation, which should not be missing in any RUMED. As part of an average concept, the use of data loggers can ensure standard-compliant monitoring and the RUMED can continue to operate. 以下產品彙編提供了一種簡單的批次控制可能性，任何 RUMED 都不應缺少這些。作為一般概念的一部分，使用數據記錄器可以確保符合標準的監控，並且 RUMED 可以繼續運行。

Included in the set 包含在套裝中:

- 1x EBI 12-TP237 壓力記錄器
- 1x AL 101 矽膠保護盒
- 1x IF 150 啟始介面
- 1x Winlog.med 軟體
- 1x PHX 800 基礎版pH測試儀(聲音信號)
- 1x AL 128 Systainer case 工具箱
- 1x AL 1100 Foam insert 泡棉鑲嵌
- 1x TDS 3 基礎版電導度測試儀

Technical data

- Measurement range 量測範圍: 0°C ...+140°C 1mbar... 4.000mbar
- Accuracy 精準度: ± 0,3°C ± 20mbar
- Resolution 解析度: 0,1°C 1mbar
- Battery 電池: Lithium 3,6V, user exchangeable 鋰 3,6V，用戶可更換




Washer disinfectors - Validation

清洗消毒器 - 驗證

Performance Qualification 性能認證 PQ

Rinsing off all process chemicals / Verification by μS measurement

沖洗掉所有流程化學品/通過 μS 測量進行驗證



The washer-disinfector manufacturer needs to be informed by the manufacturer of the chemicals about all requirements for a safe handling, data of maximum permissible residual concentrations on the products and the proof procedure used for determining the process residues. 化學品製造商需要告知清洗消毒器製造商關於安全處理的所有要求、產品上最大允許殘留濃度的數據以及用於確定流程殘留物的證明程序。

It is important to measure the conductivity.
測量電導度很重要。

Washer disinfectors - Validation

清洗消毒器 - 驗證

Rinsing off all process chemicals / Verification by μ S measurement
沖洗掉所有流程化學品/通過 μ S 測量進行驗證



EBI 12-TC230 壓力記錄器plus外加
AL 132 Flow adapter 流量適配器



Washer disinfectors - Validation

清洗消毒器 - 驗證

Disinfection performance 消毒性能

The A₀ concept 概念

The A₀-value of a disinfection procedure with moist heat is the mortification, indicated as time equivalent in seconds at a temperature of 80°C transferred to the product during the procedure, in relation to microorganisms with a Z-value of 10. 濕熱消毒程序的 A₀ 值是以 Z 值為 10 的微生物相關的在 80°C 溫度下轉移到產品的時間，以秒為單位表示

The formula for the A-value calculation is A 值計算公式為:

$$A_0 = \sum 10^{\frac{(T-80)}{Z}} \cdot \Delta t$$




Washer disinfectors - Validation

清洗消毒器 - 驗證

Thermal disinfection according to EN ISO 15883, 6.8.2

根據 EN ISO 15883, 6.8.2 進行熱消毒



1...close to the sensor for automatic control
靠近傳感器進行自動控制

2...place where the necessary temperature is reached the fastest
最快達到所需溫度的地方

3...place where the necessary temperature is reached the slowest
達到所需溫度最慢的地方

4,5,6...reference probes for the chamber temperature
腔室溫度參考探頭

It is recommended to check 2 cycles with 6 sensors or 3 cycles with 4 sensors for every loading. **The use of biological indicators instead of dataloggers is not acceptable. (ISO15883)** 建議每次裝載檢查 2 個循環（6 個傳感器）或 3 個循環（4 個傳感器）。**使用生物指示劑代替數據記錄器是不可接受的。 (ISO15883)**



Washer disinfectors - Validation

清洗消毒器 - 驗證

A_0 -values for a series of time-temperature relations
一系列時間-溫度關係的 A_0 值:

Holding time 保持時間		Temperature 溫度	A_0 -value 值
min	s	°C	
1	-	80	60
-	6	90	60
10	-	80	600
100	-	70	600
1	-	90	600
1	-	93	1 200



Most European countries use A_0 value of 3000.
大多數歐洲國家使用 3000 的 A_0 值。



How often should people perform routine monitoring?

人們應該多久進行一次例行監測？

The frequency of the routine control is the responsibility of the operator 例行控制的頻率由操作員負責。

Washer-disinfectors with temperature documentation:

Actual value with separate sensors /

control-independent recording 帶有溫度記錄的清洗消毒器：帶有單獨傳感器的實際值 / 獨立於控制的記錄

annually 每年

Washer-disinfectors with temperature documentation:

Actual value with separate sensors /

control-dependent recording 帶有溫度記錄的清洗消毒器：帶有單獨傳感器的實際值 / 依賴於控制的記錄

every 6 months 每6個月

Washer-disinfectors with temperature documentation:

Actual value without separate sensors 帶有溫度記錄的清洗消毒器：沒有單獨傳感器的實際值

every 3 months 每3個月

Washer-disinfectors without temperature documentation or with nominal value display only 沒有溫度記錄或僅顯示公稱值的清洗消毒器

daily 每日

Washer disinfectors - Validation

清洗消毒器 - 驗證

Suitable temperature and pressure logger 合適的溫度和壓力記錄器

EBI 12 temperature logger
 with 1.5 mm flex probes
 radial and axial

EBI 12 溫度記錄器
 具備 1.5 mm 柔性探頭
 徑向和軸向



EBI 11 temperature logger with
 1.5 mm flex probes EBI 11 溫度
 記錄器具備 1.5 mm 柔性探頭



EBI 12 pressure logger with 2 mm rigid probe with
 luerlock (diameter 1.95 mm)

EBI 12 壓力記錄器具備魯爾鎖 2 mm 剛性探頭（直
 徑 1.95 mm）

EBI 12 pressure logger with 1.5
 mm flex probe with luerlock

EBI 12 壓力記錄器具備 1.5 mm
 柔性探頭，具備魯爾鎖



EBI 11 temperature logger
 with 2 mm rigid probes
 (diameter 1.95 mm)


EBI 11 溫度記錄器具備 帶
 2 mm 剛性探頭
 （直徑 1.95 毫米）



EBI 11 pressure logger with
 luerlock and without luerlock
 EBI 11 壓力記錄器具備或不具
 備魯爾鎖

Washer disinfectors - Validation

清洗消毒器 - 驗證



5. Success Stories

成功的案例

Validation References -- Manufacturers and Importers

驗證參考-製造商和進口商

Worldwide contracts全球合約:

MELAG – for last two decades在過去的二十年裡

GETINGE

Miele - since自從 2002



Success in Asia在亞洲取得成功:

Sirona Dentsply, Australia New Zealand澳洲 紐西蘭

BELIMED

GETINGE Healthcare醫療保健, Australia澳洲

Miele

MMM, various countries各國



Big Hospitals – In-house validation engineers 大醫院—內部驗證工程師 SL 3111 and EBI 16



Catharina Ziekenhuis, NL 荷蘭 醫院

Covering an area of 800 square metres 佔地面積800平方米,

- 5 large chamber disinfectors 5 個大腔室消毒器
- 1 large chamber system for cleaning and disinfection 1 個用於清潔和消毒的大腔室系統
- 5 large sterilizers from Miele Professional 5 個來自德國美諾 Miele Professional 的大型消毒器



Koç University hospital, TR 土耳其 大學醫院

Covering an area of 1,000 square metres 佔地面積1000平方米,

- 12 large chamber disinfectors
12個大腔消毒器,
- 12 large sterilizers from GETINGE 12 台 GETINGE 大型消毒器



6. Solution for Other Devices

其他設備解決 方案

Plasma sterilizer process H₂O₂

電漿態過氧化氫滅菌流程

Legal background法律背景


- **ISO/DIS 22441:2021-07**
 - Sterilization of health care products – General requirements for characterization of vaporized hydrogen peroxide at low temperature 保健產品的滅菌— 低溫汽化過氧化氫特性的一般要求
 - Validation and routine control of a sterilization process for medical devices 醫療器械滅菌過程的驗證和常規控制



Plasma sterilizer process H₂O₂ 電漿態過氧化氫滅菌流程

Example of a full plasma cycle 完整的電漿態循環示例

- | | |
|-------------------|--------------------|
| 1. Vacuum 真空 | 6. Diffusion 擴散 |
| 2. Pre-plasma 前電漿 | 7. Vacuum 真空 |
| 3. Diffusion 擴散 | 8. Post-plasma 後電漿 |
| 4. Vacuum 真空 | 9. Ventilation 通風 |
| 5. Plasma 電漿 | |



Suitable temperature / pressure logger for measuring in Plasma sterilizer H₂O₂

在電漿態過氧化氫滅菌流程適合測量的溫度/壓力記錄器

EBI 12-TPX9X



描述

- Data logger for gas sterilization processes

氣體滅菌過程的數據記錄器

- The EBI 12-TP X9X pressure / temperature logger is developed for a precise measurement of pressure and temperature in low pressure sterilization processes. monitoring in H₂O₂ (plasma), formaldehyde sterilization processes

EBI 12-TP X9X 壓力/溫度記錄儀專為精確測量低壓滅菌過程中的壓力和溫度而開發。
H₂O₂ (電漿態過氧化氫)、甲醛滅菌過程中的監測

技術數據

- Measurement range 測量範圍 0°C ...+85°C
0,1mbar...1050mbar
- Accuracy精確度:
 - ± 0,1°C
 - ± 0,2mbar (0,1...50mbar)
 - ± 5% of measured value量測值 (50...100mbar) 1% FS (remaining range剩餘範圍)

Advantage from the ebro data loggers and validation equipment

ebro 數據記錄器和驗證設備的優勢

- ✓ Data logger available for process validation according to ISO 數據記錄器可用於根據 ISO 進行的流程驗證
- ✓ More than 20 years of experience with validation equipment 超過 20 年的驗證設備經驗
- ✓ TÜV certified validation equipment TÜV 認證的驗證設備
- ✓ FDA 21 CFR Part 11 / IQ,OQ,PQ documents available 可提供 FDA 21 CFR Part 11 / IQ、OQ、PQ 文件
- ✓ KompetenzCentrum in Ingolstadt 位在因戈施塔特的能力中心
- ✓ Member of the WFHSS(World Federation For Hospital Sterilisation Sciences 世界醫院消毒科學聯合會)的會員
- ✓ Member of the hygiene magazine aseptica衛生雜誌無菌的會員
- ✓ Member of the DIN(Deutsches Institut für Normung德國標準化研究所)的會員



I WILL PRACTISE my profession with conscience and dignity and in accordance with good medical practice. 我將憑良心和尊嚴並按照良好的醫療實踐從事我的職業。

WMA DECLARATION OF GENEVA, 2018
WMA 日內瓦宣言, 2018

Thank you 謝謝您!

Welcome to send questions to
歡迎發送問題至:



11471 台北市內湖區新明路273巷6號1樓
1F., No.6, Ln. 273, Xinming Rd., Neihu Dist., Taipei City 11471, Taiwan (R.O.C.)

服務專線Tel : (02)8792-3722

服務傳真Fax : (02)8792-3761

電子信箱Email : info@grandever-biotech.com.tw

公司網址Website: www.grandever-biotech.com.tw



Note: 此文件中文翻譯的部分若有進一步疑問，請參考原文或洽詢大久生物科技。

-ebro-
a xylem brand